

Super Clean™ Gas Filters

Simply the best


LC-MS System

Special LC/MS N2 Filter System
(Also available with moisture cartridges)





To meet the high flow needs of the LC/MS system, the hydrocarbon or moisture filled cartridges are positioned and placed in parallel. The incoming gas stream is split equally between the cartridges and the two streams are rejoined after purification but before the gas exits the base plate.

Without SGT Filter:



LC-MS: Sulfacetamide analysis

With SGT Filter:



www.sgt-nl.com/wheretobuy
To get more information
on where to buy our LC/MS Filters.